



# IptC Antibody

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA360016XA01SZB  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | P0ADW2   |
| Immunogen           | Recombinant Shigella flexneri IptC protein                                       |
| Raised In           | Rabbit   |
| Species Reactivity  | Shigella flexneri  |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                     |
| Form                | Liquid   |
| Conjugate           | Non-conjugated   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified  |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Product Type        | Polyclonal Antibody  |
| Immunogen Species   | Shigella flexneri  |
| Lead Time           | Made-to-order (14-16 weeks)  |
| Gene Names          | IptC   |