



luxS Antibody

Product Code	CSB-PA693377XA01FUS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Z5X1
Immunogen	Recombinant Vibrio campbellii (strain ATCC BAA-1116 / BB120) luxS protein
Raised In	Rabbit
Species Reactivity	Vibrio campbellii (strain ATCC BAA-1116 / BB120)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vibrio campbellii (strain ATCC BAA-1116 / BB120)
Lead Time	Made-to-order (14-16 weeks)
Target Names	luxS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.