



Ivr Antibody

Product Code	CSB-PA862753XA01FUT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9LBG2
Immunogen	Recombinant Leifsonia aquatica (Corynebacterium aquaticum) Ivr protein
Raised In	Rabbit
Species Reactivity	Leifsonia aquatica (Corynebacterium aquaticum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Leifsonia aquatica (Corynebacterium aquaticum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Ivr
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.