



PRXQ Antibody

Product Code	CSB-PA878722XA01FUU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9MB35
Immunogen	Recombinant Sedum lineare (Needle stonecrop) PRXQ protein
Raised In	Rabbit
Species Reactivity	Sedum lineare (Needle stonecrop)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sedum lineare (Needle stonecrop)
Lead Time	Made-to-order (14-16 weeks)
Target Names	PRXQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.