



pdxH Antibody

Product Code	CSB-PA916846XA01MXA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21159
Immunogen	Recombinant Myxococcus xanthus pdxH protein
Raised In	Rabbit
Species Reactivity	Myxococcus xanthus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Myxococcus xanthus
Lead Time	Made-to-order (14-16 weeks)
Target Names	pdxH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.