



SERPINF1 Antibody

Product Code	CSB-PA021084XA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36955
Immunogen	Recombinant Homo sapiens (Human) SERPINF1 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	SERPINF1