




otsB Antibody

Product Code	CSB-PA329093XA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31678
Immunogen	Recombinant Escherichia coli(strain K12) Trehalose-6-phosphate phosphatase (1-266aa)
Species Reactivity	Escherichia coli(strain K12)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli(strain K12)
Target Names	otsB

Image



Western Blot

Positive WB detected in Recombinant protein


All lanes: otsB antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44.5 kDa

Observed band size: 45 kDa



Western Blot

Positive WB detected in: Escherichia coli lysate

All lanes: OTSB antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.