



ospA Antibody

Product Code	CSB-PA363424XA01FUJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A3N6
Immunogen	Recombinant Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi) ospA protein
Raised In	Rabbit
Species Reactivity	Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ospA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.