



ospC Antibody

Product Code	CSB-PA600068XA01FUJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DW54
Immunogen	Recombinant Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi) ospC protein
Raised In	Rabbit
Species Reactivity	Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ospC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.