



oxIT Antibody

Product Code	CSB-PA681418XA01OFY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q51330
Immunogen	Recombinant Oxalobacter formigenes oxIT protein
Raised In	Rabbit
Species Reactivity	Oxalobacter formigenes
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oxalobacter formigenes
Lead Time	Made-to-order (14-16 weeks)
Target Names	oxIT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.