



Basic phospholipase A2 Antibody

Product Code	CSB-PA366151XA01GGN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04417
Immunogen	Recombinant Gloydius blomhoffii (Mamushi) (Agkistrodon halys blomhoffii) Basic phospholipase A2 protein
Raised In	Rabbit
Species Reactivity	Gloydius blomhoffii (Mamushi) (Agkistrodon halys blomhoffii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gloydius blomhoffii (Mamushi) (Agkistrodon halys blomhoffii)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.