



pcpB Antibody

Product Code	CSB-PA337122XA01SJK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42535
Immunogen	Recombinant Sphingobium chlorophenicum pcpB protein
Raised In	Rabbit
Species Reactivity	Sphingobium chlorophenicum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sphingobium chlorophenicum
Lead Time	Made-to-order (14-16 weeks)
Target Names	pcpB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.