



cadA Antibody

Product Code	CSB-PA890179XA01FVA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9RQB4
Immunogen	Recombinant Helicobacter felis cadA protein
Raised In	Rabbit
Species Reactivity	Helicobacter felis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Helicobacter felis
Lead Time	Made-to-order (14-16 weeks)
Gene Names	cadA