




hoc Antibody

Product Code	CSB-PA323826XA01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18056
Immunogen	Recombinant Enterobacteria phage T4 (Bacteriophage T4) hoc protein (137-376aa)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4 (Bacteriophage T4)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Target Names	hoc

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: hoc Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 34 kDa

Observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.