



hyuP Antibody

Product Code	CSB-PA156429XA01FVC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D6R8X8
Immunogen	Recombinant Microbacterium liquefaciens (Aureobacterium liquefaciens) hyuP protein
Raised In	Rabbit
Species Reactivity	Microbacterium liquefaciens (Aureobacterium liquefaciens)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Microbacterium liquefaciens (Aureobacterium liquefaciens)
Lead Time	Made-to-order (14-16 weeks)
Target Names	hyuP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.