



hyuE Antibody

Product Code	CSB-PA312426XA01PVR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00924
Immunogen	Recombinant Pseudomonas sp. (strain NS671) hyuE protein
Raised In	Rabbit
Species Reactivity	Pseudomonas sp. (strain NS671)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas sp. (strain NS671)
Lead Time	Made-to-order (14-16 weeks)
Target Names	hyuE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.