



IFNB1 Antibody

Product Code	CSB-PA011048XA01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05012
Immunogen	Recombinant Equus caballus (Horse) IFNB1 protein
Raised In	Rabbit
Species Reactivity	Equus caballus (Horse)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Equus caballus (Horse)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IFNB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.