



iaaH Antibody

Product Code	CSB-PA524365XA01EKP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O50173
Immunogen	Recombinant Enterobacter agglomerans (Erwinia herbicola) (Pantoea agglomerans) iaaH protein
Raised In	Rabbit
Species Reactivity	Enterobacter agglomerans (Erwinia herbicola) (Pantoea agglomerans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacter agglomerans (Erwinia herbicola) (Pantoea agglomerans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	iaaH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.