



hly Antibody

Product Code	CSB-PA338400XA01AXG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46028
Immunogen	Recombinant Actinobacillus pleuropneumoniae (Haemophilus pleuropneumoniae) hly protein
Raised In	Rabbit
Species Reactivity	Actinobacillus pleuropneumoniae (Haemophilus pleuropneumoniae)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Actinobacillus pleuropneumoniae (Haemophilus pleuropneumoniae)
Lead Time	Made-to-order (14-16 weeks)
Target Names	hly
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.