



hopP Antibody

Product Code	CSB-PA239353XA01FVI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A9ZM27
Immunogen	Recombinant Halomonas sp. hopP protein
Raised In	Rabbit
Species Reactivity	Halomonas sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Halomonas sp.
Lead Time	Made-to-order (14-16 weeks)
Gene Names	hopP