



HSS1 Antibody

Product Code	CSB-PA881915XA01FVJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9M4B0
Immunogen	Recombinant Senecio vulgaris (Common groundsel) HSS1 protein
Raised In	Rabbit
Species Reactivity	Senecio vulgaris (Common groundsel)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Senecio vulgaris (Common groundsel)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	HSS1