



IFNG Antibody

Product Code	CSB-PA686314XA01FVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5I6S9
Immunogen	Recombinant Macaca thibetana (Pere David's macaque) (Tibetan macaque) IFNG protein
Raised In	Rabbit
Species Reactivity	Macaca thibetana (Pere David's macaque) (Tibetan macaque)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Macaca thibetana (Pere David's macaque) (Tibetan macaque)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IFNG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.