



ibpA Antibody

Product Code	CSB-PA475840XA01ENW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B6I3S0
Immunogen	Recombinant Escherichia coli (strain SE11) ibpA protein
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain SE11)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain SE11)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ibpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.