



PCBC Antibody

Product Code	CSB-PA361478XA01FVL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05189
Immunogen	Recombinant Acremonium chrysogenum (Cephalosporium acremonium) PCBC protein
Raised In	Rabbit
Species Reactivity	Acremonium chrysogenum (Cephalosporium acremonium)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Acremonium chrysogenum (Cephalosporium acremonium)
Lead Time	Made-to-order (14-16 weeks)
Target Names	PCBC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.