



IF2CP Antibody

Product Code	CSB-PA348196XA01EUN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57997
Immunogen	Recombinant Phaseolus vulgaris (Kidney bean) (French bean) IF2CP protein
Raised In	Rabbit
Species Reactivity	Phaseolus vulgaris (Kidney bean) (French bean)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Phaseolus vulgaris (Kidney bean) (French bean)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IF2CP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.