



IFI44 Antibody

Product Code	CSB-PA011014XA01EQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27473
Immunogen	Recombinant Pan troglodytes (Chimpanzee) IFI44 protein
Raised In	Rabbit
Species Reactivity	Pan troglodytes (Chimpanzee)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pan troglodytes (Chimpanzee)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IFI44
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.