



exol Antibody

Product Code	CSB-PA291048XA01DOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E1ABX2
Immunogen	Recombinant Aspergillus ficuum exol protein
Raised In	Rabbit
Species Reactivity	Aspergillus ficuum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aspergillus ficuum
Lead Time	Made-to-order (14-16 weeks)
Gene Names	exol