



pcbC Antibody

Product Code	CSB-PA318239XA01SOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12438
Immunogen	Recombinant Streptomyces microflavus (Streptomyces lipmanii) pcbC protein
Raised In	Rabbit
Species Reactivity	Streptomyces microflavus (Streptomyces lipmanii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces microflavus (Streptomyces lipmanii)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	pcbC