



iscS Antibody

Product Code	CSB-PA015769XA01ENP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1IWD1
Immunogen	Recombinant Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks) iscS protein
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Lead Time	Made-to-order (14-16 weeks)
Target Names	iscS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.