



kamA Antibody

Product Code	CSB-PA894012XA01FVO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9XBQ8
Immunogen	Recombinant Clostridium subterminale kamA protein
Raised In	Rabbit
Species Reactivity	Clostridium subterminale
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium subterminale
Lead Time	Made-to-order (14-16 weeks)
Gene Names	kamA