



kbaZ Antibody

Product Code	CSB-PA467062XA01ENW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B6I1K5
Immunogen	Recombinant Escherichia coli (strain SE11) kbaZ protein
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain SE11)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain SE11)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	kbaZ