



waaA Antibody

Product Code	CSB-PA674948XA01DRZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q46222
Immunogen	Recombinant Chlamydia pneumoniae (Chlamydophila pneumoniae) waaA protein
Raised In	Rabbit
Species Reactivity	Chlamydia pneumoniae (Chlamydophila pneumoniae)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Chlamydia pneumoniae (Chlamydophila pneumoniae)
Lead Time	Made-to-order (14-16 weeks)
Target Names	waaA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.