



pyk Antibody

Product Code	CSB-PA339407XA01LMM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34038
Immunogen	Recombinant Lactobacillus delbrueckii subsp. bulgaricus pyk protein
Raised In	Rabbit
Species Reactivity	Lactobacillus delbrueckii subsp. bulgaricus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lactobacillus delbrueckii subsp. bulgaricus
Lead Time	Made-to-order (14-16 weeks)
Target Names	pyk
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.