



pyk Antibody

Product Code	CSB-PA343132XA01FVP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51182
Immunogen	Recombinant Sporosarcina psychrophila (Bacillus psychrophilus) pyk protein
Raised In	Rabbit
Species Reactivity	Sporosarcina psychrophila (Bacillus psychrophilus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sporosarcina psychrophila (Bacillus psychrophilus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pyk
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.