




kdsB Antibody

Product Code	CSB-PA471458XA01EOE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B5YT52
Immunogen	Recombinant Escherichia coli O157:H7 (strain EC4115 / EHEC) kdsB protein (1-248aa)
Species Reactivity	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Research Area	Others
Target Names	kdsB

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: kdsB antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.