



# kdpB Antibody

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA487115XA01EOG                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | B7N9U0                                                                           |
| Immunogen           | Recombinant Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC) kdpB protein    |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)                             |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                     |
| Form                | Liquid                                                                           |
| Conjugate           | Non-conjugated                                                                   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified                                                        |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Product Type        | Polyclonal Antibody                                                              |
| Immunogen Species   | Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)                             |
| Lead Time           | Made-to-order (14-16 weeks)                                                      |
| Gene Names          | kdpB                                                                             |