



waaA Antibody

Product Code	CSB-PA436888XA01BIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q45374
Immunogen	Recombinant Bordetella pertussis waaA protein
Raised In	Rabbit
Species Reactivity	Bordetella pertussis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bordetella pertussis
Lead Time	Made-to-order (14-16 weeks)
Target Names	waaA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.