



kefB Antibody

Product Code	CSB-PA693383XA01FRT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0ZAH7
Immunogen	Recombinant Alkalimonas amylolytica kefB protein
Raised In	Rabbit
Species Reactivity	Alkalimonas amylolytica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Alkalimonas amylolytica
Lead Time	Made-to-order (14-16 weeks)
Gene Names	kefB