



IcdC Antibody

Product Code	CSB-PA198611XA01CMC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G3KIM5
Immunogen	Recombinant Clostridium propionicum IcdC protein
Raised In	Rabbit
Species Reactivity	Clostridium propionicum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium propionicum
Lead Time	Made-to-order (14-16 weeks)
Target Names	IcdC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.