



IcnB Antibody

Product Code	CSB-PA330456XA01LNF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35518
Immunogen	Recombinant Lactococcus lactis subsp. cremoris (Streptococcus cremoris) IcnB protein
Raised In	Rabbit
Species Reactivity	Lactococcus lactis subsp. cremoris (Streptococcus cremoris)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lactococcus lactis subsp. cremoris (Streptococcus cremoris)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IcnB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.