



LBII Antibody

Product Code	CSB-PA674860XA01VFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q43296
Immunogen	Recombinant Vigna unguiculata (Cowpea) LBII protein
Raised In	Rabbit
Species Reactivity	Vigna unguiculata (Cowpea)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vigna unguiculata (Cowpea)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	LBII