



LALBA Antibody

Product Code	CSB-PA012724XA01GO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00712
Immunogen	Recombinant Capra hircus (Goat) LALBA protein
Raised In	Rabbit
Species Reactivity	Capra hircus (Goat)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Capra hircus (Goat)
Lead Time	Made-to-order (14-16 weeks)
Target Names	LALBA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.