



lat Antibody

Product Code	CSB-PA843707XA01FSR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01767
Immunogen	Recombinant Streptomyces clavuligerus (strain ATCC 27064 / DSM 738 / JCM 4710 / NBRC 13307 / NCIMB 12785 / NRRL 3585 / VKM Ac-602) lat protein
Raised In	Rabbit
Species Reactivity	Streptomyces clavuligerus (strain ATCC 27064 / DSM 738 / JCM 4710 / NBRC 13307 / NCIMB 12785 / NRRL 3585 / VKM Ac-602)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces clavuligerus (strain ATCC 27064 / DSM 738 / JCM 4710 / NBRC 13307 / NCIMB 12785 / NRRL 3585 / VKM Ac-602)
Lead Time	Made-to-order (14-16 weeks)
Target Names	lat
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.