



mlrA Antibody

Product Code	CSB-PA520643XA01SXC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E1WH94
Immunogen	Recombinant Salmonella typhimurium (strain SL1344) mlrA protein
Raised In	Rabbit
Species Reactivity	Salmonella typhimurium (strain SL1344)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Salmonella typhimurium (strain SL1344)
Lead Time	Made-to-order (14-16 weeks)
Target Names	mlrA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.