



menB Antibody

Product Code	CSB-PA708213XA01FLB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HH38
Immunogen	Recombinant Staphylococcus aureus (strain COL) menB protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus (strain COL)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Staphylococcus aureus (strain COL)
Lead Time	Made-to-order (14-16 weeks)
Target Names	menB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.