



mhpB Antibody

Product Code	CSB-PA325492XA01CUS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17295
Immunogen	Recombinant Cupriavidus necator (Alcaligenes eutrophus) (Ralstonia eutropha) mhpB protein
Raised In	Rabbit
Species Reactivity	Cupriavidus necator (Alcaligenes eutrophus) (Ralstonia eutropha)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Cupriavidus necator (Alcaligenes eutrophus) (Ralstonia eutropha)
Lead Time	Made-to-order (14-16 weeks)
Target Names	mhpB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.