



mdcC Antibody

Product Code	CSB-PA521869XA01KBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O32712
Immunogen	Recombinant Klebsiella pneumoniae mdcC protein
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Klebsiella pneumoniae
Lead Time	Made-to-order (14-16 weeks)
Target Names	mdcC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.