



mleA Antibody

Product Code	CSB-PA669680XA01OAO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q48796
Immunogen	Recombinant Oenococcus oeni (Leuconostoc oenos) mleA protein
Raised In	Rabbit
Species Reactivity	Oenococcus oeni (Leuconostoc oenos)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oenococcus oeni (Leuconostoc oenos)
Lead Time	Made-to-order (14-16 weeks)
Target Names	mleA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.