



mrpC Antibody

Product Code	CSB-PA890145XA01BRB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9RGZ3
Immunogen	Recombinant Bacillus pseudofirmus (strain OF4) mrpC protein
Raised In	Rabbit
Species Reactivity	Bacillus pseudofirmus (strain OF4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus pseudofirmus (strain OF4)
Lead Time	Made-to-order (14-16 weeks)
Target Names	mrpC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.