



mrpD Antibody

Product Code	CSB-PA870865XA01BRB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9RGZ2
Immunogen	Recombinant Bacillus pseudofirmus (strain OF4) mrpD protein
Raised In	Rabbit
Species Reactivity	Bacillus pseudofirmus (strain OF4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus pseudofirmus (strain OF4)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	mrpD